

Alarm history

The Alarm History screen displays a chronological log of all received notifications and alarms within the application. It provides a quick reference for reviewing past events and their details.




Figure 8: Alarm History screen




Figure 8: Alarm History with options menu open

Header Area

Element

Title

Options Menu (?)

Description

"Alarm History" displayed at the top of the screen. Three vertical dots in the top-right corner open an overflow menu.

Options Menu

Action

Mark All Read

Description

Marks every active alarm in the list as read, clearing the blue unread highlight on all cards.

Event List Area

The central part of the screen contains alarm cards ordered chronologically, with the most recent event at the top.

Card Element

Group name

Description

Bolded text indicating the alarm source (e.g., *Sauna*

Card Element	Description
Description	<i>Call, Nurse Call</i> . A brief description of the alarm event.
Timestamp	Date and time the event occurred, shown in the bottom-right of the card.
Unread indicator	Blue card background denotes an unread alarm.

User Flow

Tapping on any alarm card performs two actions simultaneously:

1. Opens the **Alarm Details** screen with full information about the event.
2. Marks that specific message as **Read** in the list.